



**MODERN**  
Dental Laboratory USA

# Modern PFM

## Porcelain Fused to Metal Crowns

Your team at Modern Dental brings quality, value, and customer care to everything we do – and that means PFM crowns that fit predictably, look beautiful, and have exceptional strength.

By combining the most advanced materials with technical skill and traditional craftsmanship, Modern PFM crowns provide the beauty and function you expect and your patients deserve.



*We proudly offer a  
**5-Year Warranty**  
on all Modern PFMs.*



# Modern PFM

## Porcelain Fused to Metal Crowns

### Modern PFM Advantages

Each Modern PFM substructure is carefully waxed and designed for excellent strength, utilizing alloys that are certified as biocompatible.

The metal substructures are then handstacked by our highly skilled ceramists and oven fired with premium brand ceramics. The result is optimal fit, excellent strength, and natural esthetics.



- Biocompatible
- No nickel or beryllium
- Indicated for single crowns to full arch restorations

#### High Noble

Our high noble alloys are available in White or Yellow Gold. All of our frameworks are fabricated with premium alloys to ensure we satisfy your prosthetic goals.

#### Noble:

We utilize White Chrome/ Cobalt/Pd noble alloys in our Noble frameworks which provide excellent strength and durability.

#### Non-Precious:

All non-precious frameworks are fabricated with Nickel/Berilium free alloys which provides excellent tissue compatibility.



**MODERN**  
Dental Laboratory USA

[www.moderndentalUSA.com](http://www.moderndentalUSA.com)