


FULL CAST CROWN

Certificate of Authenticity

The dental restorations used for your case were constructed of FDA approved 510k materials and hand crafted in one of our global processing centers. This is your assurance of quality and safety in materials and construction.

Item	Product Name	Manufacturer	Country of Origin	FDA 510k	Description
High Noble Alloy	Academy	Ivoclar	Liechtenstein	K002872	Au 77%, Ag 12.7%, Cu 8.5%, <1% Pt, In, Ir, Zn, Ta
High Noble Alloy	XL-X	Ivoclar	Liechtenstein	K061851	Au 62.8%, Ag 16.1%, Cu 15.7%, Pd 3.9%, <1% Ir, Zn
Noble Alloy	Minigold	Ivoclar	Liechtenstein	K905326	Au 40.0%, Ag 47.0%, Cu 7.5%, Pd 4.0 %, Zn 1.0%, <1% In, Ir
Non-Precious Alloy	Wirobond C	Bego	Germany	K021136	Co 61.0%, Cr 26.0%, Mo 6.0%, W 5.0%, <2% Si, Fe, Ce

**** This product is made entirely from one of the above alloys ****



MODERN
Dental Laboratory USA

Ramzy Abdallah, CDT, MDT
Director of Business Development