


The dental restorations used for your case were constructed of FDA approved 510k materials and hand crafted in our North American Digital Processing Center. This is your assurance of quality and safety in materials and construction.

Item	Product Name	Manufacturer	Country of Origin	FDA 510k	Description
Alloy, Base Metal	Wirobond C	Bego	Germany	K021136	Co 61%, Cr 26%, Mo 6%, W 5%, <2% Si, Fe, Ce
Porcelain	GC Initial MC, AL and LF	GC America	United States	K040567	Ceramic Powder made from SiO ₂ , Al ₂ O ₃ , K ₂ O, Na ₂ O, CaO, BaO, B ₂ O ₃ , <1% LiO ₂ , ZrO ₂ , other oxides (CeO ₂ , CeF ₃)



MODERN
Dental Laboratory USA

Ramzy Abdallah, CDT, MDT
Director of Business Development